

INTRODUCTION

Methylphenidate (MPH) is a common prescription medication that works to combat attention-deficit hyperactivity disorder (ADHD), but may also be abused recreationally. Though MPH has two chiral centers, *d*-threo-MPH is responsible for the pharmaceutical effects. Few studies have analyzed methylphenidate and its metabolites, ritalinic acid (RA) and ethylphenidate (EPH), in blood. Stability studies are crucial in a forensic setting to provide insight on ideal storage conditions and analysis time for forensic samples. This is the first method to our knowledge that analyzes the long-term stability of *d,l*-MPH, *d,l*-EPH and RA in blood.

MATERIALS & METHODS

Sample	• 250 μ L whole blood • 25 μ L internal standard solution
Preparation	• 1 mL 100 mM phosphate buffer (pH 6) • Centrifuge (2000 rpm); 10 min
Condition SPE	• 1 mL methanol • 1 mL 100 mM phosphate buffer (pH 6)
Load	• UCT Clean Screen DAU Columns (3 mL)
Wash	• 1 mL acetic acid (0.1 M) • 1 mL methanol
Elute	• 2 mL ammonium hydroxide (2%) in methanol (v/v)
Reconstitute	• Evaporate under nitrogen • Reconstitute in 100 μ L of 2:98 mobile phase

Instrumentation

- Agilent 1290 Infinity II Liquid Chromatograph
- Agilent 6470 Triple Quadrupole MS

Column

- Agilent Poroshell Chiral-V (2.1 x 100 mm, 2.7 μ m)

Mobile phase

- A: diH₂O
- B: 0.025% Ammonium acetate (w/v) + 0.0125% trifluoroacetic acid (v/v) in methanol

Flow rate

- 0.6 mL/min (isocratic, 2% A: 98% B)

RESULTS




Figure 1. MPH, EPH and RA stability at the LQC (15 ng/mL) at room temperature (A), refrigerated temperature (B), frozen temperature (C) and elevated temperature (D).




Figure 3. MPH, EPH and RA stability at the HQC (150 ng/mL) at room temperature (A), refrigerated temperature (B), frozen temperature (C) and elevated temperature (D).

Table 1. Summary of data for MPH, EPH and RA stability at the LQC (15 ng/mL) and HQC (150 ng/mL).

	Analyte Change % Loss/Increase (Timepoint deemed unstable) ¹							
	Room temp. (~25°C)		Refrigerated (4°C)		Frozen (-20°C)		Elevated (35°C)	
	LQC	HQC	LQC	HQC	LQC	HQC	LQC	HQC
<i>d</i> -MPH	-18.1 (24h)	-20.1 (24h)	-50.5 (2wk)	-54.8 (2wk)	-16.7 (5mo)	-14.3 (5mo)	-66.1 (24h)	-64.1 (24h)
<i>l</i> -MPH	-20.6 (24h)	-41.5 (24h)	-50.9 (2wk)	-56.0 (2wk)	-16.4 (5mo)	-16.9 (5mo)	-70.7 (24h)	-68.7 (24h)
<i>d</i> -EPH	-22.3 (48h)	-21.0 (48h)	-41.4 (2wk)	-47.1 (2wk)	-13.6 (5mo)	-13.4 (5mo)	-26.7 (24h)	-27.2 (24h)
<i>l</i> -EPH	-28.8 (72h)	-35.4 (72h)	-34.8 (2wk)	-53.5 (2wk)	-18.5 (5mo)	-12.1 (5mo)	-47.5 (48h)	-46.8 (48h)
RA	+52.8 (24h)	+35.8 (24h)	+40.8 (1wk)	+36.6 (1wk)	+18.5 (5mo)	+17.7 (5mo)	+116.7 (24h)	+143.5 (24h)

¹Data are displayed as %difference from 100% with the corresponding timepoint at which the analyte was deemed unstable.




Figure 2. MPH and RA concentration vs time at the LQC (15 ng/mL MPH) at room temperature (A), refrigerated temperature (B), frozen temperature (C) and elevated temperature (D).




Figure 4. MPH and RA concentration vs time at the HQC (150 ng/mL MPH) at room temperature (A), refrigerated temperature (B), frozen temperature (C) and elevated temperature (D).

Table 2. Summary of data for MPH only at the LQC (15 ng/mL) and HQC (150 ng/mL).

	Analyte Change % Loss/Increase (Timepoint deemed unstable) ²							
	Room temp. (~25°C)		Refrigerated (4°C)		Frozen (-20°C)		Elevated (35°C)	
	LQC	HQC	LQC	HQC	LQC	HQC	LQC	HQC
<i>d</i> -MPH	-25.9 (48h)	-30.4 (48h)	-26.9 (2wk)	-28.7 (1wk)	+6.5 (1mo)	+1.3 (1mo)	-53.1 (24h)	-61.1 (24h)
<i>l</i> -MPH	-34.5 (48h)	-32.6 (48h)	-29.6 (1wk)	-29.6 (1wk)	+2.8 (6mo)	+0.8 (1mo)	-60.3 (24h)	-67.4 (24h)
RA	+7.2 (1wk)	+15 (24h)	-	+6.2 (48h)	-	-	+7.2 (1wk)	+15 (24h)

²For MPH, data are displayed as %difference from 100% with the corresponding timepoint at which the analyte was deemed unstable. For RA, data are displayed as a %difference from 0% due to no RA being present in the T₀ samples. The timepoint at which RA became quantifiable is indicated.

MATERIALS & METHODS

Blood with sodium fluoride and potassium oxalate preservatives was fortified with *d,l*-MPH, *d,l*-EPH and RA at 15ng/mL and 150ng/mL. Analytes quantified with validated method.

Storage Conditions:

- Room temperature (25°C)
- Refrigerated (4°C)
- Frozen (-20°C)
- Elevated (35°C)

Time points (n=3):

- 24h, 48h, 72h, 1wk, 2wk, 3wk, 4wk, 6wk, 2mo, 3mo, 4mo, 5mo, 6mo, 7mo, 8mo, 9mo

CONCLUSIONS

- Results indicate that -20°C is the optimal storage temperature for these analytes in blood as they remained stable for up to 9 months.
- Methylphenidate degrades to ritalinic acid under non-frozen conditions.
- This study displays the importance of stability studies, gives useful information on storage conditions, and proves that quantitative values may be inaccurate due to breakdown products.

ACKNOWLEDGEMENTS

The authors would like to thank the Department of Forensic Science at Sam Houston State University for providing the resources necessary to conduct this research.

REFERENCE

Smith CR and Swortwood MJ. Short- and Long-Term Stability of Methylphenidate and Its metabolites in Blood. *Journal of Analytical Toxicology* (2021). <https://doi.org/10.1093/jat/bkab063>.