

Digital Walkthrough of the Validation of a Method for Measuring Magnetic Flux of Toner-Printed Documents

Carrie Polston*, BA, Patrick Buzzini, PhD

*Department of Forensic Science
Sam Houston State University
Huntsville, TX, USA*




Sam Houston
State University

Disclaimer

- This research was not commercially sponsored, and the authors have no financial affiliations with any commercial product manufacturers


Introduction

- Magnetic flux is a measure, in Webers, of the density of magnetic field lines passing through the surface of an object
- Magnetism has been assessed as part of the standard forensic examination scheme for toners since 1987
 - Initially proposed by the FBI as a categorical classifier
 - Updates in sensor technology allowed for the quantitation of magnetic flux



Introduction

- **Toner printed documents are increasing in prevalence in FDE casework**
 - Szafarska et al (2011) found up to 59% of a modern caseload consists of printed documents
- **Mass device production and conserved use of the electrophotographic process complicate analysis**
 - Reduce sources of variation in samples
- **Current methods provide a limited ability to discriminate**
 - Current guideline recommends 'Complementary Analyses'
 - Costly instrumentation and reagents
 - Time consuming analysis and interpretation
 - Destructive



Introduction

- Previous Research:
 - Herlaar *et al*, 2015
 - Repeatability of measurements
 - 72 samples, 19 devices
 - Biedermann *et al*, 2016
 - Repeatability and reproducibility of measurements
 - Inter-Operator variability
 - 61 samples, 61 devices
 - Mazzella and Li, 2018
 - Homogeneity of flux distribution
 - 17 samples, 6 devices and 3 cartridge types
 - Polston *et al*, 2018
 - Stability of magnetic flux fields over time
 - Relationship between toner area and flux
 - Variation in a representative population sample
 - 212 samples, 150 monocomponent toner

Introduction

- This knowledge used to facilitate method development
 - Polston *et al*, 2019
 - 5 samples, 5 devices
 - Signal suppression and enhancement
 - Hysteresis effects
 - Biasing induction current spatial effects
 - Polston *et al*, 2020
 - Reproducibility of measurements
 - 4 labs
 - Different operators, instruments, laboratory setups

Introduction


- In summary, magnetic flux of toners:
 - Is stable over time
 - Measurements are repeatable
 - Measurements can be conducted without damaging the sample
 - Data normalization is possible
 - Population variation is great enough to provide the potential to discriminate sample origin
 - Distribution is consistently inhomogeneously distributed
 - Hysteresis effects and induction spatial effects contribute significantly to variance in magnetic flux measurements
 - Must be considered in sampling methodology

Materials and Methods

- Multiple components:
 - Instrumentation
 - Validation protocol
 - Samples
 - Measurement methodology

Materials and Methods

- Instrumentation
 - Regula© Magmouse 4197
 - Used for all measurements
 - CADR Software – update version 0.98



Materials and Methods

- Validation Protocol:

- No specific QD section guidelines, so adapted relevant sections from toxicology ASB validation guidelines

Validation Parameters	Criteria
LOD	*Previously assessed - Magnetic particles visible but no measured response
LOQ	*Previously assessed - Magnetic particles visible and measurement value 1 nWb
Precision	± 20%
Bias	± 20%
Suppression/Enhancement	*Previously assessed – hysteresis effects, biasing induction current effects
Sample Stability	*Previously assessed – stable long-term
Interferences	No impact on measurements
Matrix Effects	No impact on measurements

Materials and Methods

- Samples
 - QC samples
 - Negative and positive QC samples
 - Matrix samples
 - 5 different matrix sources
 - Interference study samples
 - 25 different materials
 - Pen inks, markers, inkjet inks
 - Population sample
 - 54 toner samples
 - 20 monocomponent toners

Materials and Methods

- Sample Information

- Matrices
- Interferents
- Population

Manufacturer	Weight (g/sm)
Skilcraft	32
HP	32
Neenah	65
Georgia Pacific	32
Clairfontaine	90

Media	Number of Samples	Manufacturers	Colors	Models
Ballpoint Pens	8	6	3	7
Gel Pens	3	2	2	2
Rollerball Pens	2	2	1	2
Felt Tip Markers	4	2	3	2
Inkjet Inks	8	2	4	2

Manufacturer	Samples	Color	BW	Bicomponent	Mono
Brother	7	3	4	7	0
Canon	8	4	4	4	4
HP	27	12	15	12	15
Kyocera	1	0	1	0	1
OKI	2	0	2	2	0
Ricoh	3	1	2	3	0
Xerox	6	3	3	6	0

Materials and Methods

Measurement Methodology

- Bias and Precision
 - 2 sites on QC sample (low and high)
 - 5 replicates
 - 3 sessions
- Matrix Effects and Interferences
 - 3 sites on sample
 - 3 replicates
- Population
 - 1 site on sample
 - 5 replicates
 - ILS protocol

Results

Precision and Bias

	Toner QC data					
	Low QC			High QC		
	Day 1	Day 2	Day 3	Day 1	Day 2	Day 3
Rep 1	30	33	27	57	63	57
Rep 2	31	34	28	58	64	59
Rep 3	31	28	29	62	58	67
Rep 4	31	32	29	62	57	61
Rep 5	32	28	29	63	61	60
p-w	2.3	9.1	3.1	4.5	5.0	6.2
p-b	6.9			4.9		
	27.71	59.33				
Bias	Low	High				
	8.7	2.1				



Results

- Precision and Bias
- Interferences and Matrix Effects

Validation Parameters	Results
LOD	*Previously assessed - Magnetic particles visible but no measured response
LOQ	*Previously assessed - Magnetic particles visible and measurement value 1 nWb
Precision	2.3%-9.1%
Bias	2.1%-8.7%
Suppression/Enhancement	*Previously assessed – hysteresis effects, biasing induction current effects
Sample Stability	*Previously assessed – stable long-term
Interferences	No impact on measurements
Matrix Effects	No impact on measurements

Results

- Population



Conclusions

- Magnetic flux measurement method validation can be conducted successfully, with precision and bias within acceptable ranges
- No matrix effects or interferences were found to occur
- Magnetic flux measurements provide the potential for a robust, non-destructive toner examination method

References

- Szafarska M., Wietecha-Posluszny R., Wozniakiewicz M., Koscielniak P., Examination of color ink-jet printing inks by capillary electrophoresis, *Talanta*, 84 (2011), pp. 1234-1243.
- Herlaar, K.; Mieremet, M.; and Fakkel, M. (2015). Measuring magnetic properties to discriminate between different laser printers. *Journal of the American Society of Questioned Document Examiners*, Vol. 18.2, pp. 51-66.
- Biedermann, A.; Bozza, S.; Taroni, F.; Fürbach, M.; Li, B.; and Mazzella, W.D. (2016). Analysis and Evaluation of Magnetism of Black Toners on Documents Printed by Electrophotographic Systems. *Forensic Science International*, Vol. 267, pp. 157-165.
- Mazzella, W.D. and Li, B.B. (2018) Is Magnetic Flux a Valuable Tool for the Analysis of Electrophotographic-Printed Documents? *Journal of Forensic Science and Medicine*, Vol. 4, pp. 197-201.
- Polston C., Mazzella W., Fürbach M., Buzzini P. (2018) Assessing the Repeatability and Reproducibility of Magnetic Flux Measurements and their Potential to Discriminate Toner Printed Documents. *Journal of the American Society of Questioned Document Examiners*, Volume 21, Number 2, pp. 45-56.
- Polston C., Mazzella W., Buzzini P. "Investigating the Interlaboratory Reproducibility of Magnetic Flux Measurements of Toners." Presented at the American Academy of Forensic Sciences 72nd Annual Meeting, Anaheim, CA, February 17-22 2020.
- Polston C., Thompson D., Buzzini P. "Assessing the Impacts of Induction Spatial Effects on Magnetic Flux Measurements of Toners." Presented at the American Society of Questioned Document Examiners 77th Annual Meeting, Cary, NC, August 4-8 2019.
- Academy Standards Board ANSI/ASB Standard 036, Standard Practices for Method Validation in Forensic Toxicology (2019).

Questions

Carrie Polston

cep020@shsu.edu



Sam Houston
State University