

Exploring Immunomagnetic Isolation of Spermatozoa

Julia Wang, MS^{1*}; Samantha Davis, MS¹; Sheree Hughes, PhD¹; Brendan Chapman, PhD²; Andrew Currie, PhD²; Rachel Houston, PhD¹

¹Department of Forensic Science, Sam Houston State University, Huntsville, TX, 77340

²Medical, Molecular and Forensic Science, Murdoch University, Perth, Australia



INTRODUCTION

Forensic DNA labs often process sexual assault kits, where separating female epithelial cells from male spermatozoa is critical but time-consuming using conventional methods. These methods involve preferential lysis with reducing agents, which can lead to DNA loss¹. Biomolecule-targeting cell isolation methods, commonly used in clinical research, offer a potential solution. This study focuses on sperm-specific proteins—SPAM1, SPACA1, and ZBPB1—shown to improve sperm cell enrichment². The objective was to assess whether magnetic beads conjugated to antibodies targeting these proteins can effectively isolate sperm cells from mixtures.

MATERIALS AND METHODS

Immunocytochemistry

- 4% Formaldehyde fixation, R.T. 15 minutes
- 5% Normal goat serum blocking, R.T. 2 hours
- 1° Antibody incubation, 4°C overnight
- 2° Antibody & DAPI incubation, R.T. 1 hour
- Leica FS CB microscope

Bead Preparation & Cell Isolation

- Performed following manufacturer recommendations with Dynabeads® M-450 Epoxy (Thermo Fisher Scientific, Waltham, MA).
- Antibodies used were anti-rabbit IgG polyclonal anti-SPAM1, anti-SPACA1, and anti ZBPB (Thermo Fisher Scientific)

DNA Extraction, Quantification, & Isolation

- Bead-isolated sperm or vaginal epithelial cells
- QIAamp® DNA Investigator Kit
- Investigator® Quantiplex® Pro Kit

RESULTS

Antibody Binding Confirmation




Figure 1. Target protein binding confirmation by immunocytochemistry in spermatozoa. Each target protein was detected by anti-SPAM1, anti-SPACA1, anti-ZBPB, or combined in an antibody cocktail. Scale bars represents 50 μm.

Antibody Cocktail-Bead Specificity




Figure 2. Non-specific interaction between BSA-blocked magnetic beads and cells. Increasing incubation time increased non-specific interactions between the magnetic bead complexes and both spermatozoa and vaginal epithelial (VE) cells.

Antibody Cocktail-Bead Sensitivity




Figure 3. Minimal differences observed between control beads and antibody cocktail-conjugated beads. Majority of cells were found in supernatant fraction rather than the bead-complex fraction.

CONCLUSIONS

- Immunocytochemistry confirmed antigen-binding capability of all three antibodies (**Figure 1**).
- Vaginal epithelial cells showed cross-reactivity with the magnetic beads (**Figure 2**).
- Similar results between antibody-coated and uncoated beads indicated non-specific cell-bead interactions (**Figure 3**).

ACKNOWLEDGEMENTS

Samples were collected from donors using SHSU IRB protocol 2020-248. This work was supported by financial support of the Forensic Sciences Foundation Lucas Grant. The opinions, findings, and conclusions or recommendations expressed in this presentation are those of the authors and do not necessarily reflect those of the Forensic Sciences Foundation.

REFERENCES

1. Vuichard S, Borer U, Bottinelli M, Cossu C, Malik N, Meier V, et al. Differential DNA extraction of challenging simulated sexual-assault samples: a Swiss collaborative study. *Investig Genet* 2011;2:11. <https://doi.org/10.1186/2041-2223-2-11>.
2. Zhao X-C, Wang L, Sun J, Jiang B-W, Zhang E-L, Ye J. Isolating Sperm from Cell Mixtures Using Magnetic Beads Coupled with an Anti-PH-20 Antibody for Forensic DNA Analysis. *PLOS ONE* 2016;11(7):e0159401. <https://doi.org/10.1371/journal.pone.0159401>.